

**Q: What is a TTL®?**

**A: A Training/Test Lung.**

The TTL® :


- embodies a mechanical simulation of the human pulmonary system capable of measuring more than twenty individual parameters.
- visually demonstrates almost any normal or pathological pulmonary condition.
- provides accurate measurements of volumes, pressures and flow rates of medical equipment replacing several measuring instruments at a fraction of their combined costs.
- accommodates most gas specific measuring sensors and other sensing equipment.
- available with PneuView® 3 software or in an analog version.

**Applications:**

- Classroom Instruction
- Ventilator Testing
- Product Evaluation
- Product Demonstration
- Pulmonary Research
- Quality Control
- Design Engineering
- Clinical Intervention




**Single Adult Test Lung**




Mobility– Need a portable solution? The Single Adult lung provides compact, reliable simulation options and offers the same compliance and resistance options as our other adult lungs.

**Adult/Infant Test Lung**

Versatility– Need a wider spectrum of simulation? The Adult/Infant lung will provide an adult and infant lung in one solution. The infant lung is also ideal for HFOV and jet ventilation testing.



**Dual Adult Test Lung**



Capability– Need to measure large volumes? The Dual Adult lung can provide up to 4 liters of capacity.

Compliance for any lung can be set independently by a variable position spring located on the side. Resistance values can be selected using parabolic resistors introduced into the airway.

**Breath Simulation Module**

The Breath Simulation Module transforms the Dual Adult models into spontaneously breathing lung systems.



- Simulate a spontaneously breathing patient for training or studying the use of CPAP, IMV, SIMV and pressure support of other ventilation models designed for use with breathing patients.
- Test and troubleshoot devices designed to support a spontaneously breathing patient for training.
- Measure added work of breathing associated with breathing circuits.
- Assess trigger sensitivity and response of demand flow systems.

Breath Rate:	1-60 Breaths per Minute
Inspiratory Time:	User selectable - .50, 1.0, 1.5, & 2.0 sec (+/- .015 sec)
Breath Volume:	50-2000 mL
Manual Mode:	Manual trigger provided for individual breaths

**NOTE:** The BSM may be used on Adult/Infant TTL Models but its use is limited. It can only be used to create a spontaneously breathing infant lung.

PneuView® 3 Software

PneuView® 3 enhances the demonstration of ventilation phenomena. Possessing a wide array of capabilities and features, PneuView3® will complement any equipment testing, certification or calibration program.

Select from up to 26 data parameters gathered, calculated and reported. All parameters are available for review at any time.

- Visually demonstrates, in real-time, the relationship between pressure, volume, and flow waveforms.
- Provides acquisition, storage, and review of data.
- Trend ventilator performance for up to 1000 hours.
- Measures pressure, volume, flow and timing parameters.
- Provides FiO<sub>2</sub> and ambient temperature measurements.



PneuView® 3 Software

CALCULATIONS:

- Breath Rate
- Inspiratory Time
- Expiratory Time
- I:E Ratio
- Tidal Volume
- Minute Volume
- Baseline Pressure
- and more

Computer Requirements:

PneuView® 3 software requires minimal computer resources to ensure proper operation.

- 1 GHz processor
- 1 GB RAM (32-bit) or 2 GB RAM (64-bit)
- 1 GB available hard disk space
- DirectX 9 graphics device with WDDM 1.0 or higher driver.
- Windows® 7 (32 or 64 bit) OS or newer
- Microsoft .NET Framework 4.5 or higher

Specifications:

	Single Lung	Adult/Infant Lung		Dual Adult Lung	
	Adult	Adult	Infant	Each Lung	Total
Tidal Volume:	2.0 L	2.0 L	200 mL	2.0 L	4.0 L
Residual Lung Volume:	986 mL	986 mL	200 mL	986 mL	1.972 L
Lung Compliance: (adjustable) L/cmH <sub>2</sub> O <small>Accuracy: +/- 3% (at calibration volumes)</small>	.01 to .10	.01 to .10	.001 to .01	.01 to .10 (each lung)	NA
Airway Resistance: (selectable) cm H <sub>2</sub> O/L/sec <small>Accuracy: +/- 5% (at calibration flows)</small>	Rp5, Rp20 or Rp50	Rp5, Rp20 or Rp50	Rp50, Rp200 or Rp500	Rp5, Rp20 or Rp50	NA
Pressure: (Lung/Airway) cm H <sub>2</sub> O <small>Accuracy: +/- 0.5% of user-calibrated Full Scale Output (FSO)</small>	-20 to 120	-20 to 120	-20 to 120	-20 to 120	NA
Weight:	25 lbs	37 lbs		37 lbs	

Contact us for more information.

Phone: 616.554.9696  
Toll Free: 800.530.9939  
International: +1.616.554.9696  
E-mail: [sales@michinst.com](mailto:sales@michinst.com)

Or visit our website at [www.michiganinstruments.com](http://www.michiganinstruments.com)



Innovations Like This Don't Come Out Of Thin Air



Announcing the NEW Michigan Lung

Provides accurate simulation of a wide range of normal and diseased lung conditions for ventilator testing/calibration and respiratory therapy instruction.