



Quick Start Guide

Thank you for purchasing the PneuView®3 (PV3).

1. Insert the flash drive into any USB port.
2. Copy the PneuView3 folder onto the computer desktop.

Inside the PneuView3 folder you will find the following subfolders and their contents-

Info-promotional brochures and videos for our Training/Test lung and CPR products.

Install-setup file to install the PV3 software onto your computer.

Manuals-TTL and PV3 software operation manuals in .pdf format.

Photos-a picture of each available model to use for desktop wallpaper if desired.

Michigan Instrument .com-a shortcut to our website. Copy to your desktop for convenience.

To install the software-

Open the Install folder and double click on the setup file.

During the installation you will be asked which folder to install the software in. Also, you can choose to install a desktop shortcut (which is highly recommended). After the installation is complete, you have an option to launch the PV3 software. If you elect to do so, have the lung setup, plugged in, powered on and the USB cable plugged into the computer.